

TORS[®]

TORSIONAL ORTHOPAEDIC REVISION SYSTEM

Summit
Medical 

Summit
Medical 

TORS[®]

TORSIONAL ORTHOPAEDIC REVISION SYSTEM

DUAL MODE ULTRASOUND FOR REVISION OF ENDOPROSTHESIS

- Integrated cement removal and haemostatic soft tissue dissection
- Brand new switched handset format
- A range of new torsional probe geometries
- Performance guaranteed



TORS is distributed by:

Summit Medical Group Limited
Bourton Industrial Park
Bourton on the Water
Gloucestershire
GL54 2HQ
United Kingdom

Tel: +44 (0)1451 821311
Email: info@summit-medical.co.uk

www.summit-medical.com



RSL is the legal manufacturer of TORS
and owns all IP rights.

Radley Scientific Ltd
Bremridge House
Ashburton
Devon
TQ13 7JX

LIT239-01

www.summit-medical.com

PERFORMANCE EASE OF USE RELIABILITY

TORS is the newest ultrasonic system for bone cement removal.

The use of innovative, patented technologies enable a significantly more effective, gentler, and safer bone cement removal than comparable systems on the market.

The newly designed probes increase the effectiveness of cement removal, and thus reduces operating theatre working time.

The well-engineered system has extremely stable functionality and performance, as well as easy and comfortable handling.



CEMENT REMOVAL

The two-part TORS cement handset comprises:
Inner transducer – a simple switchless canister with robust seals.

Outer handpiece – a lightweight, ergonomic, single use plastic handle with comfy annular switch accessible at all angles.

SOFT TISSUE

The soft tissue transducer is a completely reusable device that will enable haemostatic skin and muscle incisions equivalent to monopolar diathermy but running much cooler and therefore smoke-free yet with better bleeding control.

The soft tissue transducer can be used in primaries as well.

TORS KIT



T1G
TORS Generator



T1CT
TORS Cement Transducer (Re-usable)



T1TT
TORS Soft Tissue Transducer (Re-usable)



T1FC
TORS Cement Footswitch



T1FT
TORS Soft Tissue Footswitch



T1CC
T1CC TORS Cement Cable

RE-USABLE PROBES – 200mm



P4R2
TORS Cement 200 Ø4 Piercer Probe



P6R2
TORS Cement 200 Ø6 Piercer Probe



P8R2
TORS Cement 200 Ø8 Piercer Probe



P10R2
TORS Cement 200 Ø10 Piercer Probe



S6R2
TORS Cement 200 Ø6 Scraper Probe



S8R2
TORS Cement 200 Ø8 Scraper Probe



S10R2
TORS Cement 200 Ø10 Scraper Probe

RE-USABLE PROBES – 100mm



P6R1
TORS Cement 100 Ø6 Piercer Probe



P8R1
TORS Cement 100 Ø8 Piercer Probe



S6R1
TORS Cement 100 Ø6 Scraper Probe



S8R1
TORS Cement 100 Ø8 Scraper Probe

EXTENSION BARS



ECR1
TORS Curved Extension Bar



ESR1
TORS Straight Extension Bar (Short)



ESR2
TORS Straight Extension Bar (Long)

ACCESSORIES



T1S
TORS Spanner



T1AT + T1SM
TORS Instrument Tray



T1IT
TORS Generator Cart

CONSUMABLES



T1AH
TORS Axial Handpiece



T1CW
TORS Cleaning Cell